|                                                                                                                                                    |                                                                                                                                                   |                                                                |                                                               |
|----------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|---------------------------------------------------------------|
| <br><b>PT ETANA</b><br><b>BIOTECHNOLOGIES</b><br><b>INDONESIA</b> | <b>STANDARD OPERATING PROCEDURE</b>                                                                                                               |                                                                | Page : 1 of 8                                                 |
|                                                                                                                                                    | <b>TITLE :</b><br><b>PERAWATAN BERKALA TERHADAP MESIN LAF MOBILE TYPE ORABS</b><br><b>PREVENTIVE MAINTENANCE OF LAF MOBILE MACHINE TYPE ORABS</b> |                                                                | SOP No. : SOP/EBI/EN-065<br>Revision : 00<br>Effective Date : |
| <b>PREPARED BY :</b>                                                                                                                               | <b>REVIEWED BY :</b>                                                                                                                              | <b>REVIEWED BY :</b>                                           |                                                               |
| <b>PRIBADI SETIA BAYU AJI</b><br><b>(NON UTILITY ENGINEERING TECHNICIAN)</b>                                                                       | <b>WENDI RUKMANSYAH</b><br><b>(ENGINEERING SPV)</b>                                                                                               | <b>PURNO BUDI KISWANTO</b><br><b>(ENGINEERING SR. MANAGER)</b> |                                                               |
| <b>REVIEWED BY :</b>                                                                                                                               | <b>APPROVED BY :</b>                                                                                                                              | <b>APPROVED BY :</b>                                           |                                                               |
| <b>VINO SOADUON</b><br><b>(PROD SR. MANAGER)</b>                                                                                                   | <b>DEVIKA SRI WULANDARI</b><br><b>(QS MANAGER)</b>                                                                                                | <b>SHANDY</b><br><b>(QA MANAGER)</b>                           |                                                               |

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*This SOP is used by the Engineering department as a guideline in carrying out maintenance on the LAF Mobile Type ORABS (Open Restricted Access Barrier System) machines in the production area.*

### **3. TANGGUNG JAWAB RESPONSIBILITY**

#### **3.1 Teknisi bertanggung jawab:**

*Technician are responsible:*

3.1.1. Melakukan perawatan sesuai jadwal yang telah ditetapkan.

*Perform maintenance according to a predefined schedule.*

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*Follow up Engineering Services Requests from other departments related ORABS issue.*

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*Report any deviation to Engineering Supervisor.*

#### **3.2 Supervisor Engineering bertanggung jawab:**

*Engineering Supervisor is responsible:*

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*Create schedule of maintenance periodically.*

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*Ensure implementation of maintenance is according to the schedule.*

3.2.3. Mengadakan *review* bulanan terhadap hasil kerja teknisi.

*Conduct a monthly review to works of technicians.*

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*Provide training to the related department.*

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*Coordinate and report to Engineering Manager if there is malfunction to ORABS System.*

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*Calibration & Validation supervisor is responsible:*

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*Performs validation and calibration of measuring instruments used on ORABS system like magnehelic.*

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*Engineering Manager is responsible:*

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*Determine replacement parts if the appropriate parts are no longer available for sale by proposing changes if there are changes to the parts used. Refer to point 5.7.*

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*Decide and coordinate action required, if necessary replacement, equipment modification, change system of ORABS System.*

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*Periodic maintenance carried out according to a defined schedule to prevent equipment malfunction or deterioration in performance. Activities include inspection, cleaning, calibration, and replacement of parts.*

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*The process of verifying and adjusting measuring instruments so that the measurement results comply with predetermined standards.*

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*A low-speed, uniform airflow system, parallel and perpendicular to the work surface, designed to keep the area free from airborne particle contamination.*

#### 4.4. *ORABS (Open Restricted Access Barrier System)*

*ORABS (Open Restricted Access Barrier System)*

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*An open barrier system with laminar airflow (LAF) and HEPA filters that protects critical areas from contamination.*

#### 4.5. *HEPA (High Efficiency Particulate Air filter)*

*HEPA (High Efficiency Particulate Air filter)*

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*High efficiency particulate air filter. This type of filter is suitable for use as a cleaner or for implementation in rooms with special air conditioning.*

#### 4.6. *UPS (Uninterruptible Power Supplies)*

*UPS (Uninterruptible Power Supplies)*

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*Electronic devices capable of providing temporary backup power when the main power supply is interrupted.*

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*The way of maintenance should be done regularly and in accordance with the schedule.*

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*Trained or qualified training technicians who can conduct maintenance of ORABS system.*

#### 5.3. Keselamatan kerja harus tetap di perhatikan selama perawatan.

*Work safety must be kept in mind during the maintenance.*

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*Regular maintenance activities should be documented after maintenance done.*

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*If the measuring instrument is found to be in a defective state, it should be clearly marked.*

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*Work equipment must be in accordance with the set.*

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*Spare parts replacement must be in accordance with the specifications and functions.*

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*For major underlying changes and changing the ORABS system must use a change control form and be approved by related managers.*

## 6. ALAT DAN BAHAN

### TOOLS AND MATERIAL

#### 6.1. ALAT

##### TOOLS

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*Wrench Key*

- 6.1.2. Obeng

*Screwdriver*

- 6.1.3. Kunci L

*Allen Key*

#### 6.2. BAHAN

##### MATERIAL

- 6.2.1. Ethanol 70%

*Ethanol 70%*

- 6.2.2. Kain Wipes

*Wipes fabric*


## 7. PROSEDUR

### PROCEDURE

- 7.1. Perawatan 3 bulanan (L2) pada *Open Restricted Access Barrier System (ORABS)*.

*Three (3) month maintenance 3 month (L2) of Open Restricted Access Barrier System (ORABS).*

|                                                                                                                                                    |                                                                                                                                                                           |                                               |
|----------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|
| <br><b>PT ETANA</b><br><b>BIOTECHNOLOGIES</b><br><b>INDONESIA</b> | <b>STANDARD OPERATING PROCEDURE</b>                                                                                                                                       | Page 6 of 8                                   |
|                                                                                                                                                    | <b>TITLE :</b><br><b>PERAWATAN BERKALA TERHADAP MESIN</b><br><b>LAF MOBILE TYPE OBRAS</b><br><br><b>PREVENTIVE MAINTENANCE OF LAF MOBILE</b><br><b>MACHINE TYPE OBRAS</b> | SOP No. : SOP/EBI/EN-065<br><br>Revision : 00 |

*M = Mechanical*

*E = Electrical*

Gambar 1. Mesin LAF Mobile Type ORABS (Open Restricted Access Barrier System)

Figure 1. LAF Mobile Machine Type ORABS (Open Restricted Access Barrier System)

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*Check the blower motor from abnormal sound, repair or replace if necessary. Look at the figure 1 code M01.*

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*Check the single-phase power connection and repair any if visible damage. Look at the figure 1 code E.01.*

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*Check the pressure difference of the HEPA filter. If it is less than 100 pascals, it means there is a leak in the HEPA, HEPA seal, or Magnahelic hose. If it is more than 500 pascals, it means it is clogged, and must be repaired immediately. If the HEPA pressure still does not meet the requirements (Magnahelic Pressure Requirements 200–300 pascal) after repair, immediately inform the production team to request time to replace the HEPA filter. HEPA filters should be replaced every 5 years at most if the integrity results still meet the requirements (integrity requirements < 0.005%). Look at the figure 1 code M02.*

|                                                                                                                                                    |                                                                                                                                                                           |                                               |  |
|----------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|--|
| <br><b>PT ETANA</b><br><b>BIOTECHNOLOGIES</b><br><b>INDONESIA</b> | <b>STANDARD OPERATING PROCEDURE</b>                                                                                                                                       | Page 7 of 8                                   |  |
|                                                                                                                                                    | <b>TITLE :</b><br><b>PERAWATAN BERKALA TERHADAP MESIN</b><br><b>LAF MOBILE TYPE OBRAS</b><br><br><b>PREVENTIVE MAINTENANCE OF LAF MOBILE</b><br><b>MACHINE TYPE OBRAS</b> | SOP No. : SOP/EBI/EN-065<br><br>Revision : 00 |  |

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*Check and tighten all cable terminal connections, circuit breakers, relays, replace if there is a color change on the electrical components. Look at the figure 1 code E02.*

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*Check the condition of the display controller, repair or replace it if it looks damaged. Look at figure 1 code E03.*

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*Check the wheel bolts, tighten them if they are loose and replace the wheel if they appear damaged. Look at the figure 1 code M03.*

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*Check the condition of gloves for ORABS. Glove leakage testing is performed by the production team during the glove leak test process using glove leak test (GLT). Look at the figure 1 code M04.*

- 7.1.8. Periksa kekencangan baut2 pada pengikat cover – cover di mesin.

*Check the tightness of the bolts on the cover fasteners on the machine.*

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*Check the condition of the acrylic glass. If there are any visible cracks or breaks, immediately inform the supervisor so that it can be repaired or replaced with a new one.*

## 8. REFERENSI

### REFERENCE

Petunjuk SOP/EBI/FF-027 Pengoperasian dan Pembersihan Open Restricted Access Barrier System (ORABS) Mobile.

*SOP/EBI/FF-027 Instructions for Operating and Cleaning the Mobile Open Restricted Access Barrier System (ORABS).*

|                                                                                                                                                    |                                                                                                                                                                           |                                               |
|----------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|
| <br><b>PT ETANA</b><br><b>BIOTECHNOLOGIES</b><br><b>INDONESIA</b> | <b>STANDARD OPERATING PROCEDURE</b>                                                                                                                                       | Page 8 of 8                                   |
|                                                                                                                                                    | <b>TITLE :</b><br><b>PERAWATAN BERKALA TERHADAP MESIN</b><br><b>LAF MOBILE TYPE OBRAS</b><br><br><b>PREVENTIVE MAINTENANCE OF LAF MOBILE</b><br><b>MACHINE TYPE OBRAS</b> | SOP No. : SOP/EBI/EN-065<br><br>Revision : 00 |

## 9. LAMPIRAN

### ATTACHMENT

Formulir No. SOP/EBI/EN-065-F01A Perawatan Berkala Terhadap Mesin *LAF Mobile Type ORABS*.

*Form No. SOP/EBI/EN-065-F01A Periodic Maintenance of LAF Mobile Type ORABS Machines.*

## 10. RIWAYAT

### HISTORY

| <b>Tanggal</b><br><i>Date</i> | <b>No. Change Control</b><br><i>Change Control No.</i> | <b>No. Revisi</b><br><i>No. Revision</i> | <b>Alasan Perubahan</b><br><i>Reason of Change</i> |
|-------------------------------|--------------------------------------------------------|------------------------------------------|----------------------------------------------------|
|                               | CC-EBI-EN-015-25                                       | 00                                       | Dokumen baru<br><i>New Document</i>                |

## 11. DOKUMEN TERKAIT

### RELATED DOCUMENT

Tidak ada

N/A

## 12. TARGET PELATIHAN

### TRAINING TARGET

| <b>Departemen</b><br><i>Department</i>                    | <b>Jabatan Personil</b><br><i>Personil Title</i> |
|-----------------------------------------------------------|--------------------------------------------------|
| Departement Engineering<br><i>Engineering Departement</i> | Teknisi<br><i>Technician</i>                     |

## 13. DISTRIBUSI

### DISTRIBUTION


- 13.1. Departemen Engineering (*Hard copy*).  
*Engineering Department (Hard copy).*
  - 13.2. Departemen Sistem Mutu (*Hard copy*).  
*Quality System Departement (Hard copy).*
  - 13.3. Departemen Produksi (*Hard copy*).  
*Production Department (Hard copy).*
-

|                                                                                                                                                    |                                                                                                                                                         |                                       |
|----------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------|
| <br><b>PT ETANA</b><br><b>BIOTECHNOLOGIES</b><br><b>INDONESIA</b> | <b>FORMULIR / FORM</b>                                                                                                                                  | Document No. :<br>SOP/EBI/EN-065-F01A |
|                                                                                                                                                    | <b>PERAWATAN BERKALA TERHADAP</b><br><b>MESIN LAF MOBILE TYPE ORABS</b><br><br><b>PREVENTIVE MAINTENANCE OF LAF</b><br><b>MOBILE MACHINE TYPE ORABS</b> |                                       |

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| Level | Component | Activity |   |   |   |   |   |   |   | Condition |  | Remark |
|-------|-----------|----------|---|---|---|---|---|---|---|-----------|--|--------|
|       |           | C        | I | L | R | S | T | W | B | G         |  |        |

|           |                                                     |  |  |  |  |  |  |  |  |  |  |  |
|-----------|-----------------------------------------------------|--|--|--|--|--|--|--|--|--|--|--|
| <b>L2</b> |                                                     |  |  |  |  |  |  |  |  |  |  |  |
|           | Blower motor, M01                                   |  |  |  |  |  |  |  |  |  |  |  |
|           | Connection power 1 phasa, E01                       |  |  |  |  |  |  |  |  |  |  |  |
|           | Pressure magnehelic (Requirement 200 – 300 Pa), M02 |  |  |  |  |  |  |  |  |  |  |  |
|           | Connection terminal, circuit braker and relay, E02  |  |  |  |  |  |  |  |  |  |  |  |
|           | Function Display controller, E03                    |  |  |  |  |  |  |  |  |  |  |  |
|           | Condition wheel, M03                                |  |  |  |  |  |  |  |  |  |  |  |
|           | Gloves for ORABS, M04                               |  |  |  |  |  |  |  |  |  |  |  |
|           | tightness of the bolts                              |  |  |  |  |  |  |  |  |  |  |  |
|           | condition of the acrylic glass                      |  |  |  |  |  |  |  |  |  |  |  |

|        |  |  |  |  |  |  |  |  |  |  |  |  |
|--------|--|--|--|--|--|--|--|--|--|--|--|--|
| Note : |  |  |  |  |  |  |  |  |  |  |  |  |
|--------|--|--|--|--|--|--|--|--|--|--|--|--|

| No | Date | Start | Finish | Duration |        | Total personnel | Area | Remark |
|----|------|-------|--------|----------|--------|-----------------|------|--------|
|    |      |       |        | Hour     | Minute |                 |      |        |
|    |      |       |        |          |        |                 |      |        |
|    |      |       |        |          |        |                 |      |        |

|                                                                                                                                                                                                                                                                          |                                   |                                           |                                              |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|-------------------------------------------|----------------------------------------------|
| Remark :<br>C : Cleaning (pembersihan)<br>I : Inspection (pemeriksaan)<br>L : Lubrication (pelumasan)<br>R : Replacement (penggantian)<br>S : Setting (penyetelan)<br>T : Tightened (dikencangkan)<br>W : Welding (pengelasan)<br>G : Good ( baik )<br>B : Bad ( jelek ) | Level :<br><br>L2 : Every 3 month | Done by :<br><br>(                      ) | Checked by :<br><br>(                      ) |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|-------------------------------------------|----------------------------------------------|