 PT ETANA BIOTECHNOLOGIES INDONESIA	STANDARD OPERATING PROCEDURE		Page : 1 of 11
	TITLE : PERAWATAN TERHADAP SEMUA MESIN- MESIN DI PT. ETANA BIOTECHNOLOGIES INDONESIA MAINTENANCE OF ALL MACHINES AT PT. ETANA BIOTECHNOLOGIES INDONESIA		SOP No. : SOP/EBI/EN-004 Revision : 02 Effective Date : Reviewed Date :
PREPARED BY :	PREPARED BY :	PREPARED BY :	
WENDI RUKMANSYAH (ENGINEERING SUPERVISOR)	RIKI DEPARNO (ENGINEERING SUPERVISOR)	BURHANUDDIN ROBBANI (ENGINEERING SUPERVISOR)	
REVIEWED BY :		APPROVED BY :	
PURNO BUDI KISWANTO (ENGINEERING MANAGER)		HAPPY MONDA PINTAULI (QA MANAGER)	

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This SOP is used by the Department of Engineering as a guide in handling maintenance and repairing utility, production, warehouse, quality assurance (QA), and quality control (QC) machines at PT. Etana Biotechnologies Indonesia.

3. TANGGUNG JAWAB

RESPONSIBILITY

- 3.1. Teknisi Engineering *Preventive Maintenance* bertanggung jawab :
Engineering Technician Preventive Maintenance are responsible :

- 3.1.1. Melakukan tugas perawatan yang sesuai dengan SOP ini.
Perform maintenance tasks in accordance with this SOP.
- 3.1.2. Melakukan perawatan sesuai jadwal yang telah ditetapkan.
Perform maintenance according to a predefined schedule.
- 3.1.3. Melaporkan setiap terjadi kerusakan dan perbaikan yang diperlukan kepada supervisor engineering.
Report any damage and repairs that are needed to the engineering supervisor.
- 3.2. Supervisor Engineering bertanggung jawab :
Engineering Supervisor are responsible :
- 3.2.1. Membuat jadwal perawatan secara berkala.
Make a schedule of maintenance periodically.
- 3.2.2. Memastikan pelaksanaan perawatan sesuai dengan jadwal.
Ensure the implementation of maintenance is according to the schedule.
- 3.2.3. Memastikan semua alat ukur yang terlibat dalam proses perawatan sesuai dengan spesifikasi yang telah ditentukan (dikalibrasi).
Ensure that all measuring instruments used in the maintenance process are meet the specifications (calibrated).
- 3.2.4. Mengadakan kajian bulanan terhadap hasil kerja teknisi preventive maintenance.
Perform a monthly review of the workings of the preventive maintenance technicians.
- 3.2.5. Memberikan pelatihan kepada departemen lain yang terkait.
Provide training to others department which related.
- 3.2.6. Melakukan evaluasi pada jadwal perawatan mesin tahunan
Evaluate the annual machine maintenance schedule
- 3.3. Manajer Engineering bertanggung jawab :
Engineering Manager are responsible:
- 3.3.1. Menentukan spare parts pengganti apabila spare parts yang sesuai sudah tidak dijual lagi.
Determine replace parts if suitable spare part are no longer sold.
- 3.3.2. Memutuskan dan mengkoordinir tindakan yang diperlukan, apabila perlu dilakukannya penggantian, modifikasi peralatan, perubahan sistem dan menyiapkan prosedur perubahan.
Resolve and coordinate necessary actions, if necessary replacement, equipment modification, system changes and prepared change control.
-

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Preventive maintenance is maintenance carried out on a scheduled basis, generally periodically, where a number of maintenance tasks such as inspection, repair, replacement, cleaning, lubrication and adjustments are carried out.

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Corrective maintenance is a maintenance activity carried out if or when the specified preventive maintenance schedule cannot be carried out according to schedule so it must be revised according to the agreement with the user.

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Autonomous maintenance is a machine independent maintenance performed by the machine operator. If all this time the operators are only trained to operate the machines, then it's time to be further trained. Operators should be trained to be able to detect minor irregularities in the machine and make repairs themselves.

The goal of the autonomous maintenance is to develop the operator's ability to be able to detect symptoms of damage before the actual damage occurs. For this reason the operator must first create a regular workplace so that any machine irregularities can be detected quickly

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Work safety must always be considered, use personal protective equipment (PPE) that fits its designation on conditions when handling, maintaining and repairing all machines in PT. Etana Biotechnologies Indonesia, on dangerous parts use labels or signs to avoid work errors.

- 5.2. Ceklis perawatan disusun berdasarkan buku manual, instruksi vendor, dan pengalaman di lapangan.

Maintenance checklist based on manual book, vendor instruction and experience on the field.

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The maintenance checklist on the form has several levels of maintenance, which include L1 maintenance every 1 month, L2 maintenance every 3 months, L3 maintenance every 6 months and L4 maintenance every 1 year.

For the level checklist on the maintenance form it is adjusted to the preventive schedule, if the L1 level checklist is only in the L1 column we will check it, if we check L2 level it is L1 and L2, for L3 level we check L1 L2 and L3, if it is L4 level we check it L1 L2 L3 and L4, each form also has a different inspection technique.

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Only technician who have been trained and qualified who can work carry out maintenance / modification of a production machine.

- 5.5. Jika teknisi sedang melakukan *preventive maintenance* dan terdapat hal seperti :

If technicians do preventive maintenance and there are things like :

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*If a dangerous event is found, report it to Engineering Supervisor and Engineering Manager.
If needed, report to the security officer.*

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Every spare part replacement specification must have a form change control.

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Evaluation of the annual schedule is reviewed every 1 month which is contained in the monthly schedule and if there is a new machine in the production area, warehouse and quality assurance (QC), it will be included in the schedule after qualification is carried out by revising the annual schedule. Machines that are no longer used for user needs will be excluded from the annual PM schedule.

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Preventive maintenance has a goal for :

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Reduce the possibility of damaged products.

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Periodic Maintenance Procedures in utility, production, warehouse, quality assurance (QA) and quality control (QC) machines.

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Create a Preventive Maintenance annual schedule first then ask for approval from the engineering manager, production manager, quality control (QC) manager, and quality accurance (QA) manager.

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Based on the annual preventive maintenance schedule, pour it into the monthly schedule, where there are levels (L1) once a month, (L2) once every 3 months, (L3) once every 6 months and (L4) once a year, on each machines accordance with the type of maintenance SOP that has been made and then seek approval from the engineering manager, production manager, quality control (QC) manager, and quality assurance manager (QA), where this monthly schedule can be implemented

according to schedule or change dates according to related department activities but must be carried out in the month according to the annual schedule.

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1 month before preventive maintenance is carried out, communicate / notify the relevant departments to plan their production schedule which at that time will be used for preventive maintenance, so that the implementation of preventive maintenance could be completed on time.

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Record all preventive maintenance activities in the inspection form of each machine according to the SOP on the machines, then record the summary of the implementation in the Machine History Card No. SOP/EBI/EN-004-F03.

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Supervision of utility, production, quality control (QC) and quality assurance (QA) machines during long holidays.

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For supervision process, Engineering provides form to list machines or tools that must be monitored during a long holiday and users who will fill out the list of machines or equipment that will supervised during the long holiday with the approval of the relevant supervisor and manager.

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Add a reason for the change to the history column in the document attachment no. SOP/EBI/EN-004-L01. Any changes to the annual machine maintenance schedule.

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			<i>Revision of the Annual Machine Maintenance Schedule Template with document number. : SOP/EBI/EN-004-L01B by adding point 2 about History.</i>
--	--	--	--

11. DISTRIBUSI

DISTRIBUTION

11.1. Departemen Engineering.

Engineering Department.

11.2. Departemen Produksi

Production Department.

11.3. *Department* Manajemen Material bagian Warehouse.

Warehouse Section Material Management Department

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Annual machine maintenance schedule template is made as follows:

TEMPLATE JADWAL PERAWATAN MESIN TAHUNAN
ANNUAL MACHINE MAINTENANCE SCHEDULE TEMPLATE

Document No. : SOP/EBI/EN-004-L01B

Effective Date :

[illegible]

L1	: Every 1 month
L2	: Every 3 month
L3	: Every 6 month
L4	: Every 1 year

Tanggal Date	No. Change Control Change Control No.	No. Revisi No. Revision	Alasan Perubahan Reason of Change

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Monthly machine maintenance schedule template is made as follows:

TEMPLATE JADWAL PERAWATAN MESIN BULANAN
MONTHLY MACHINE MAINTENANCE SCHEDULE TEMPLATE

Document No. : SOP/EBI/EN-004-L02B

Effective Date :

Department: Engineering

[illegible]

Remark :

Level Service :

Prepared by :

Approved by :

L1	Every 1 month
L2	Every 3 month
L3	Every 6 month
L4	Every 1 year

Technician

Eng. Supervisor

Val & Cal. Supervisor

Eng. Manager

Prod. Manager

PPIC. Manager

QA. Manager

 PT ETANA BIOTECHNOLOGIES INDONESIA	FORMULIR / FORM		Document No. : SOP/EBI/EN-004-F03C
	KARTU RIWAYAT MESIN MACHINE HISTORY CARD		Effective Date :

Machine Name :

Location :

Date	Working Code	Machine Parts	No. PJE	Repair information	Spare parts are replaced	Performance status			Done by	Checked by
						Running Time (Min)				
						10	30	60		

Working code information :

A : Break Down Maint 1 = Repair
 B : Preventive Maint 2 = Replacement
 3 = Improvement