 PT ETANA BIOTECHNOLOGIES INDONESIA	STANDARD OPERATING PROCEDURE		Page : 1 of 8	
	TITLE : PERAWATAN TERHADAP MESIN DEBAGGING TOFFLON MAINTENANCE OF DEBAGGING MACHINE TOFFLON		SOP No. : SOP/EBI/EN-048 Revision : 00 Effective Date : Reviewed Date :	
PREPARED BY :	REVIEWED BY :	APPROVED BY :	APPROVED BY :	
PRIBADI SETIA BAYU AJI (NON UTILITY ENGINEERING TECHNICIAN)	INDRA PRASTOWO (NON UTILITY ENGINEERING SPV)	PURNO BUDI KISWANTO (ENGINEERING MANAGER)	HAPPY MONDA PINTAULI (QA MANAGER)	

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Provides a guidance for maintenance of the Debagging machine Tofflon in accordance with standards at PT. Etana Biotechnologies Indonesia.

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This SOP is used by Engineering Department especially for technician as guidance in order to maintain of Debagging machine Tofflon in the production 1st floor.

3. TANGGUNG JAWAB

RESPONSIBILITY

3.1. Teknisi Engineering *Preventive Maintenance* bertanggung jawab :

Engineering Technician Preventive Maintenance are responsible :

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Performs maintenance activities properly, in accordance with this SOP.

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Performs maintenance of Debagging machine Tofflon according to a predefined schedule.

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Report any damage and repairs that are needed to the Engineering Supervisor.

3.2. Supervisor Engineering bertanggung jawab :

Engineering Supervisor are responsible :

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Make sure the implementation of maintenance is according to the schedule.

3.2.3. Mengadakan *review* bulanan terhadap hasil kerja Teknisi *Preventive Maintenance*.

Perform a monthly review of the work results of the preventive maintenance technicians.

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Provide training to related personnel.

3.3. Manajer Engineering bertanggung jawab :

Engineering Manager are responsible :

3.3.1 Menentukan *spare part* pengganti apabila *spare part* yang sesuai tidak dijual lagi.

Determining replacement of spare part if the appropriate spare part are no longer sold.

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Resolve and coordinate necessary actions, if necessary replacement, equipment modification or system changes and prepared change control document.

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PM (Preventive Maintenance) is scheduled routine maintenance to ensure equipment or assets continue to function normally according to their function and to prevent unexpected damage in the future

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The air filters can filter particles up 99,995% with size 0,30 µm.

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6.1.4 Tangga lipat

Ladder

Commented [i3]: English nya Multimeter?

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Check all machine equipment to operate properly

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Check the electrical circuit and sensor. Make sure the power is the same

With the standart and the sensor can read the nest.

7.1.1.3. Pengecekan pisau pemotong. Pastikan pisau tajam

Checking the cutting system. Make sure the knife is sharp.

7.1.1.4. Pemeriksaan koneksi selang pneumatic

Checking the connection of pneumatic tube.

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Check the cutter surface, repair or changing the cutter if needed. Caution when removing the cutter, don't touch sharp edge.

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Changing the pre filter if needed, for changing the pre filter follow this procedure :

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Prepare ladder, plastic bag and wrench 8.

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Open the bolt on the cover of pre filter, do it carefully

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11. DISTRIBUSI

DISTRIBUTION

11.1. Departemen Pemastian Mutu
Quality Assurance Department

11.2. Departemen Engineering
Engineering Department

11.3. Departemen Produksi
Production Department

11.4. Departemen Pengendali Kualitas
Quality Control Department

 PT ETANA BIOTECHNOLOGIES INDONESIA	FORMULIR / FORM PERAWATAN BERKALA TERHADAP MESIN DEBAGGING TOFFLON	Document No. : SOP/EBI/EN-048-F01A
	PERIODIC MAINTENANCE OF DEBAGGING MACHINE TOFFLON	Effective Date :

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Level	Component	Activity							Condition		Remark
		C	I	L	R	S	T	W	B	G	
L2											
	Pisau Pemotong/Cutter										
	Electrical current and sensor										
	Pneumatic tubing connection										
L4											
	Check the HEPA Filter										
	Check the Pre Filter										
	Electrical Connection										

Note :

No	Date	Start	Finish	Duration		Total personnel	Area	Remark
				Hour	Minute			

Remark : C : Cleaning (pembersihan) I : Inspection (pemeriksaan) L : Lubrication (pelumasan) R : Replacement (penggantian) S : Setting (penyetelan) T : Tightened (dikencangkan) W : Welding (pengelasan) G : Good (baik) B : Bad (jelek)	Level : L2 : Every 3 month L4 : Yearly	Done by : ()	Checked by : ()
--	--	---	--