 PT ETANA BIOTECHNOLOGIES INDONESIA	STANDARD OPERATING PROCEDURE	Page : 1 of 6
	TITLE : PERAWATAN TERHADAP MESIN WATER BATH DAIHAN SCIENTIFIC MAINTENANCE OF WATER BATH DAIHAN SCIENTIFIC MACHINE	SOP No. : SOP/EBI/EN-050 Revision : 00 Effective Date : Reviewed Date :
PREPARED BY :	REVIEWED BY :	APPROVED BY :
 WENDI RUKMANSYAH (NON UTILITY ENGINEERING SPV)	 PURNO BUDI KISWANTO (ENGINEERING MANAGER)	 HAPPY MONDA PINTAULI (QA MANAGER)

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Provide clear instructions on how to maintenance for the Water Bath Daihan Scientific Machine in the production area.

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This SOP is used by the Engineering department as a guideline for carrying out Machine Waterbath Daihan Scientific maintenance in the production area.

3. TANGGUNG JAWAB

RESPONSIBILITY

3.1. Teknisi Engineering *Preventive Maintenance* bertanggung jawab :

Engineering Technician Preventive Maintenance are responsible :

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Performing maintenance duties properly, in accordance with this SOP.

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Perform maintenance according to a predefined schedule.

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Report any damage and repairs that are needed to the Engineering Supervisor.

3.2. Supervisor Engineering bertanggung jawab :

Engineering Supervisor are responsible :

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the machine remains good, cleaned every machine that is the responsibility, especially on the inside of the crust that is attached.

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Ensure the condition of electrical and electronic equipment that is part of the machine placed in water-prone areas, in order to always avoid splashing water.

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Conduct condensate water discharge at Air Service Regulator on machine equipped with pneumatic system.

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The Calibration and Validation Staff from Department Quality Assurance is responsible for calibrating the components of the measuring instrument used, eg Thermometer, Pressure Gauge, Thermo Controller etc, and qualify regularly.

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The Daihan Scientific Water Bath Machine is a laboratory tool in the form of a container used to heat liquids by immersing them in water.

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Autonomous maintenance is an independent engine maintenance carried out by the machine operator. If all this time the operators are only trained to operate the engine, then it's time to be further trained.

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Before implementing preventive maintenance, engineering staff must coordinate with the operator about the area to prepare the area.

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When carrying out maintenance, the related operator must accompany Engineering to avoid operating errors when the machine is running or the operator can undergo autonomous maintenance on the machine.

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Work safety must be paid attention to at the time of preventive maintenance, use PPE (Personal protective equipment) which according to allotment in condition when maintenance.

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For the replacement of critical components or spare parts that have a direct impact on the product when replaced with different brands and different types but one type, make change control.

6. ALAT DAN BAHAN

TOOLS AND MATERIALS

6.1. ALAT

TOOLS

1 set tools box

1 set of toolbox

6.2. BAHAN

MATERIALS

6.2.1. Ethanol 20%.

Ethanol 20%.

6.2.2. Kain Wipes..

Wipes Fabric.

7. PROSEDUR

PROCEDURE

Perawatan 3 bulanan (L2) pada Mesin Water Bath Daihan Scientific.

Maintenance 1 month (L2) of Water Bath Daihan Scientific Machine.

Gambar 1. Mesin Water Bath Daihan Scientific

Figure 1. Waterbath Daihan Scientific Machine

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Check the connection of the ON / OFF switch on the machine, if there is a loose one, fix it or replace it immediately.

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Check the potential connection on the machine, if something is loose, fix it or replace it immediately.

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Check the condition of the other switch connections, if something is loose, fix it or replace it immediately.

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Check the power cable connections, tighten if any are loose, and replace the power cable if one appears to be damaged.

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Clean the machine body if there is dirty.

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Check the inside of the machine, open the cover, and observe, if there are things that have the potential to cause damage, do repairs immediately.

8. REFERENSI

REFERENCE

Manual book of periodic maintenance Waterbatch Daihan Scientific Machine.

9. LAMPIRAN

ATTACHMENT

Formulir perawatan berkala terhadap Mesin Waterbatch Daihan Scientific List No. SOP/EBI/EN-050-F01A.

Periodically Maintenance form of Machine Waterbatch Daihan Scientific List No. SOP/EBI/EN-050-F01A.

10. RIWAYAT

HISTORY

Tanggal <i>Date</i>	No. Change Control <i>Change Control No.</i>	No. Revisi <i>No. Revision</i>	Alasan Perubahan <i>Reason of Change</i>
	CC-EBI-EN-011-23	00	Dokumen baru <i>New Document</i>

11. DISTRIBUSI

DISTRIBUTION

11.1. Departemen Engineering.

Engineering Department.

11.2. Departemen Pemastian Mutu (QA).

Quality Assurance Department.

11.3. Departemen Produksi.

Production Department.

 PT ETANA BIOTECHNOLOGIES INDONESIA	FORMULIR / FORM PERAWATAN BERKALA PADA MESIN WATER BATH DAIHAN SCIENTIFIC	Document No. : SOP/EBI/EN-050-F01A
	PERIODIC MAINTENANCE OF MACHINE WATER BATH DAIHAN SCIENTIFIC	Effective Date :

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Level	Component	Activity							Condition		Remark
		C	I	L	R	S	T	W	B	G	
L2											
	Connection of the ON / OFF switch										
	Potential connection										
	Condition of the other switch connections										
	Power cable connections										
	Body machine										
	Check the inside of the machine										

Note :

No	Date	Start	Finish	Duration		Total personnel	Area	Remark
				Hour	Minute			

Remark : C : Cleaning (pembersihan) I : Inspection (pemeriksaan) L : Lubrication (pelumasan) R : Replacement (penggantian) S : Setting (penyetelan) T : Tightened (dikencangkan) W : Welding (pengelasan) G : Good (baik) B : Bad (jelek) * : Critical component	Level : L2 : Every 3 month	Done by : ()	Checked by : ()
--	-----------------------------------	---	--