 PT ETANA BIOTECHNOLOGIES INDONESIA	STANDARD OPERATING PROCEDURE	Page : 1 of 7
	TITLE : PERAWATAN TERHADAP MESIN LEAK TEST LT-PRO HA 3570 BONFIGLIOLI MAINTENANCE OF LEAK TEST LT- PRO HA 3570 BONFIGLIOLI MACHINE	SOP No. : SOP/EBI/EN-054 Revision : 00 Effective Date : Reviewed Date :
PREPARED BY :	REVIEWED BY :	APPROVED BY :
 WENDI RUKMANSYAH (NON UTILITY ENGINEERING SPV)	 PURNO BUDI KISWANTO (ENGINEERING MANAGER)	 HAPPY MONDA PINTAULI (QA MANAGER)

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Provide clear instructions on how to maintenance for the Leak Test LT-PRO HA 3570 Bonfiglioli Machine in the production area.

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This SOP is used by the Engineering department as a guideline for carrying out Machine Leak Test LT-PRO HA 3570 Bonfiglioli maintenance in the production area.

3. TANGGUNG JAWAB

RESPONSIBILITY

3.1. Teknisi Engineering *Preventive Maintenance* bertanggung jawab :

Engineering Technician Preventive Maintenance are responsible :

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Performing maintenance duties properly, in accordance with this SOP.

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Perform maintenance according to a predefined schedule.

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Report any damage and repairs that are needed to the Engineering Supervisor.

3.2. Supervisor Engineering bertanggung jawab :

Engineering Supervisor are responsible :

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the machine remains good, cleaned every machine that is the responsibility, especially on the inside of the crust that is attached.

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Ensure the condition of electrical and electronic equipment that is part of the machine placed in water-prone areas, in order to always avoid splashing water.

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Conduct condensate water discharge at Air Service Regulator on machine equipped with pneumatic system.

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The Calibration and Validation Staff from Department Quality Assurance is responsible for calibrating the components of the measuring instrument used, eg Thermometer, Pressure Gauge, Thermo Controller etc, and qualify regularly.

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LT-PRO HA 3570 Leak Testing Machine is a laboratory instrument used to check for leaks in the primary packaging of a product by using pressure or vacuum within a certain range of values.

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Autonomous maintenance is an independent engine maintenance carried out by the machine operator. If all this time the operators are only trained to operate the engine, then it's time to be further trained.

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Before implementing preventive maintenance, engineering staff must coordinate with the operator about the area to prepare the area.

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When carrying out maintenance, the related operator must accompany Engineering to avoid operating errors when the machine is running or the operator can undergo autonomous maintenance on the machine.

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Work safety must be paid attention to at the time of preventive maintenance, use PPE (Personal protective equipment) which according to allotment in condition when maintenance.

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For the replacement of critical components or spare parts that have a direct impact on the product when replaced with different brands and different types but one type, make change control.

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
Check the sound or vibration of the vacuum when the engine is running, then observe, if the sound of the vacuum is abnormally harsh, replace the vacuum motor bearing.

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Check the performance of the vacuum machine, the standard parameters can be seen in the production SOP.

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Check the oring at the top for the connector, if the oring looks damaged or the condition is not normal, immediately replace the oring.

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Check the air filter section on the engine body, do the cleaning, if the filter is cleaned it still looks dirty replace the filter.

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Check the oring on the other connecting parts, if the oring looks damaged or the condition is not normal, immediately replace the oring.

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Check the electrical connection parts, fix if there is a loose one.

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Check the sensor for the pressure, clean it, if necessary we recalibrate it.

8. REFERENSI

REFERENCE

Manual book of periodic maintenance Leak Test LT-PRO HA 3570BONFIGLIOLI

9. LAMPIRAN

ATTACHMENT

Formulir perawatan berkala terhadap Mesin Leak Test LT-PRO HA 3570 Bonfiglioli List No. SOP/EBI/EN-054-F01A.

Periodically Maintenance form of Machine Leak Test LT-PRO HA 3570 Bonfiglioli List No. SOP/EBI/EN-054-F01A.

10. RIWAYAT

HISTORY

Tanggal <i>Date</i>	No. Change Control <i>Change Control No.</i>	No. Revisi <i>No. Revision</i>	Alasan Perubahan <i>Reason of Change</i>
	CC-EBI-EN-014-23	00	Dokumen baru <i>New Document</i>

11. DISTRIBUSI

DISTRIBUTION

11.1. Departemen Engineering.

Engineering Department.

11.2. Departemen Pemastian Mutu (QA).

Quality Assurance Department.

11.3. Departemen Produksi.

Production Department.

 PT ETANA BIOTECHNOLOGIES INDONESIA	FORMULIR / FORM PERAWATAN BERKALA PADA MESIN LEAK TEST LT-PRO HA 3570 BONFIGLIOLI	Document No. : SOP/EBI/EN-054-F01A
	PERIODIC MAINTENANCE OF MACHINE LEAK TEST LT-PRO HA 3570 BONFIGLIOLI	Effective Date :

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Level	Component	Activity							Condition		Remark
		C	I	L	R	S	T	W	B	G	
L2											
	Vibration of the vacuum										
	Performance of the vacuum										
	Check the oring at the top										
	Air filter										
	Oring on the other connecting										
	Electrical connection										
	Sensor for the pressure										

Note :

No	Date	Start	Finish	Duration		Total personnel	Area	Remark
				Hour	Minute			

Remark : C : Cleaning (pembersihan) I : Inspection (pemeriksaan) L : Lubrication (pelumasan) R : Replacement (penggantian) S : Setting (penyetelan) T : Tightened (dikencangkan) W : Welding (pengelasan) G : Good (baik) B : Bad (jelek) * : Critical component	Level : L2 : Every 3 month	Done by : ()	Checked by : ()
--	-----------------------------------	---	--