 PT ETANA BIOTECHNOLOGIES INDONESIA	STANDARD OPERATING PROCEDURE		Page : 1 of 12	
	TITLE : PERAWATAN TERHADAP MESIN AUTOMATIC VISUAL INSPECTION TRUKING MAINTENANCE AUTOMATIC VISUAL INSPECTION TRUKING OF MACHINE		SOP No. : SOP/EBI/EN-061 Revision : 00 Effective Date : Reviewed Date :	
PREPARED BY :	REVIEWED BY :	APPROVED BY :	APPROVED BY :	
WENDI RUKMANSYAH (ENGINEERING SPV)	PURNO BUDI KISWANTO (ENGINEERING MANAGER)	VINO SOADUON (PROD MANAGER)	SHANDY (QA MANAGER)	

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Provide clear instructions on how to maintenance for Truking Automatic Visual Inspection Machine in the production area.

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This SOP is used by the Engineering department as a guideline for carrying out Truking Automatic Visual Inspection Machine maintenance in the production area.

3. TANGGUNG JAWAB

RESPONSIBILITY

3.1. Teknisi Engineering *Preventive Maintenance* bertanggung jawab :

Engineering Technician Preventive Maintenance are responsible :

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Performing maintenance duties properly, in accordance with this SOP.

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Perform maintenance according to a predefined schedule.

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Report any damage and repairs that are needed to the Engineering Supervisor.

3.2. Supervisor Engineering bertanggung jawab :

Engineering Supervisor are responsible :

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Make schedule of maintenance periodically.

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Ensure implementation of maintenance is according to the schedule.

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Ensure all measuring instruments involved in the maintenance process are in accordance with the specifications specified (calibrated).

3.2.4. Mengadakan pengkajian bulanan terhadap hasil kerja Teknisi Preventive Maintenance.

Perform a monthly review of the workings of the Preventive Maintenance Technicians.

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Provide training to others department which related.

3.3. Manajer Engineering bertanggung jawab :

Engineering Manager is responsible :

3.3.1. Menentukan spare parts pengganti apabila spare parts yang sesuai sudah tidak dijual lagi.

Determine replace parts if suitable spare part are no longer sold.

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Resolve and coordinate necessary actions, if necessary replacement, equipment modification, system changes and prepared change control.

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The Production Machine Operator is responsible for :

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Report if the area is detected by defects / abnormalities that normally impede the operation of the machine.

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Perform daily maintenance such as: tighten ferulle clamps loose bolts, clean the body of the machine/Main Frame at all times so that it does not rust and the surface of the machine remains good, cleaned every machine that is the responsibility, especially on the inside of the crust that is attached.

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Ensure the condition of electrical and electronic equipment that is part of the machine placed in water-prone areas, in order to always avoid splashing water.

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Conduct condensate water discharge at Air Service Regulator on machine equipped with pneumatic system.

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Calibration and Validation Staff from Department Quality Assurance is responsible for calibrating the components of the measuring instrument components used, eg Thermometer, Pressure Gauge, Thermo Controller etc, and perform qualification regularly.

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Automatic Visual Inspection Machine is an automatic system that uses visual technology such as cameras and optical sensors to carry out automatic inspections of products or components, which aims to detect defects or nonconformities in products quickly and accurately during the production process or after production is complete.

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Autonomous maintenance is one of the pillars of Total Productive Maintenance (TPM), which aims to empower machine operators to perform basic maintenance independently, without always relying on maintenance technicians, which aims to Increase operator involvement in maintaining equipment and Reduce downtime and improve machine efficiency.

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Preventive maintenance is a type of maintenance that is carried out on a scheduled and planned basis to prevent damage or machine failure before it occurs. The main objective is to keep the equipment in optimal condition, avoid sudden breakdowns, and extend the life of the machine.

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Before implementing preventive maintenance, engineering staff must coordinate with the operator about the area to prepare the area.

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When carrying out maintenance, the related operator must accompany Engineering to avoid operating errors when the machine is running or the operator can undergo autonomous maintenance on the machine.

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Work safety must be paid attention to at the time of preventive maintenance, use PPE (Personal protective equipment) which according to allotment in condition when maintenance. Such as safety helmets, safety shoes, and gloves if the maintenance position is outside the production building.

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For the replacement of critical components or spare parts that have a direct impact on the product when replaced with different brands and different types but one type, make change control.

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Staff or operators who perform preventive maintenance must be trained and understand how to perform maintenance correctly.

6. ALAT DAN BAHAN

TOOLS AND MATERIALS

6.1. ALAT

TOOLS

- 6.1.1. 1 set *tools box*, pada saat melakukan *Preventive Maintenance*.

1 set of toolbox, when performing Preventive Maintenance

- 6.1.2. APD (Alat Pelindung Diri) seperti *sepatu safety, helm safety dan sarung tangan safety*

PPE (Personal Protective Equipment) such as safety shoes, safety helmets and safety gloves

6.2. BAHAN

MATERIAL

- 6.2.1. Kain Bebas Serat

Lint Free Fabric

- 6.2.3. Pelumas *Spray Food Grade*

Food Grade Lubricant spray

7. PROSEDUR

PROCEDURE

Perawatan 3 bulanan (L2) dan 1 Tahunan (L4) pada Mesin *Automatic Visual Inspection* Truking.
Maintenance 3 month (L2) and 1 Year (L4) of Automatic Visual Inspection Truking Machine.

- 7.1. Perawatan 3 bulanan (L2).

3 month Maintenance (L2).

Gambar 1. Bagian electrical

Figure 1. Electrical elements

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Maintenance of Electrical Elements.

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Check the servo drive battery, if there is a failure or alarm of other electrical elements during operation, immediately replace the battery, battery replacement is done every 3 years.

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Check electrical power connections and other connections.

7.1.2. Perawatan pada unit vacuum pump becker.

Maintenance on the becker vacuum pump unit

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
Clean the body unit vacuum and vacuum pump room using fabric material. Look at the figure 2 code M.01

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Clean the air filter in the unit vacuum using compressed air.

7.1.3. Perawatan pada *gears* dibagian bawah mesin. *Look at the figure 2 code M.02*
Maintenance of the gears at the bottom of the machine.

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
Clean with a lint-free fabric and regreasing the new one.

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Give lubricant to the clamp claws, using food grade spray lubricant.

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Press on the clamp claws, if there are any that are hard, repair them so they run smoothly again.

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Maintenance of section pneumatic system.

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Clean the condensation water form compressed air using lint-free fabric.

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
Clean the air filter of section regulator compressed air using lint-free fabric.

7.1.5.3. Bersihkan solenoid–solenoid pada komponen *pneumatic* menggunakan kain bebas serat.

Clean the solenoid of pneumatic component using lint-free fabric.

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Check, if dirty change section air filter of tubing hose pneumatic.

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Lubrication of moving parts such as hinges, joints and bearings to reduce friction and wear. Use food grade lubricant.

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Clean the robotic arm from dust, dirt and debris that can interfere with operations, avoid using harsh chemicals that can damage components.

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Check the condition of cables and connectors to ensure there is no damage, wear, or loose connections. Repair or replace damaged cables immediately.

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Check the condition of the motor and gearbox to ensure there is no excessive wear or damage, if there is any indication of damage, immediately repair or replace it.

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Check for unusual sounds or vibrations that could indicate a problem, and correct the vibrations.

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
Check the rubber parts of the suction cup, replace them if any are torn or damaged.

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Check the starwheel fastener, make sure it is tight so that the starwheel remains balanced with the other starwheels.

Gambar 8. Vacuum pump becker

Figure 8. Becker vacuum pump

7.2. Perawatan 1 Tahunan (L4).

Maintenance 1 year (L4)

7.2.1. Perawatan pada unit Vacuum pump becker.

Maintenance on the Becker vacuum pump unit.

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Check and clean the carbon vanes, replace if damaged. (The position of carbon vanes is inside the vacuum, need to make a cover to check carbon vanes). Look at the figure 8 code M.03.

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Check section rubber coupling, replace if there broken. Look at the figure 8 code M.03.

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Check air filter, replace if there broken. Look at the figure code M.02

8. REFERENSI

REFERENCE

Buku manual Perawatan berkala Truking Mesin Inspeksi Visual Otomatis.

Manual book of periodic Maintenance Truking Automatic Visual Inspection Machine.

9. LAMPIRAN**ATTACHMENT**

Formulir perawatan berkala terhadap Mesin *Automatic Visual Inspection Truking* SOP/EBI/EN-061-F01A.

Periodically Maintenance form of Truking Automatic Visual Inspection Machine List No. SOP/EBI/EN-061-F01A.

10. RIWAYAT**HISTORY**

Tanggal <i>Date</i>	No. Change Control <i>Change Control No.</i>	No. Revisi <i>No. Revision</i>	Alasan Perubahan <i>Reason of Change</i>
	CC-EBI-EN-014-24	00	Dokumen baru <i>New Document</i>

11. DISTRIBUSI**DISTRIBUTION**

11.1. Departemen Engineering.

Engineering Department.

11.2. Departemen Pemastian Mutu (QA).

Quality Assurance Department.


11.3. Departemen Produksi.

Production Departement.


11.4. Departemen Sistem Mutu.

Quality System Departement.

 PT ETANA BIOTECHNOLOGIES INDONESIA	FORMULIR / FORM PERAWATAN TERHADAP MESIN AUTOMATIC VISUAL INSPECTION TRUKING	Document No. : SOP/EBI/EN-061-F01A
	MAINTENANCE OF AUTOMATIC VISUAL INSPECTION TRUKING OF MACHINE	Effective Date :

Level :											
Month :											
Date :											
Day :											
Hours Meter :											
Level	Component	Activity							Condition		Remark
		C	I	L	R	S	T	W	B	G	
L2	Electrical										
											
Servo drive											
Connection power and other											
Vacuum pump (Becker)											
											
Body vacuum											
Air filter											
Gears (Bottom of the machine)											
											
Cleaning and regreasing											
Lubricant for AS gear											
Clamp claws											
											
Lubricant to the clamp claws											
Press on the clamp claws											

 PT ETANA BIOTECHNOLOGIES INDONESIA	FORMULIR / FORM PERAWATAN TERHADAP MESIN AUTOMATIC VISUAL INSPECTION TRUKING MAINTENANCE OF AUTOMATIC VISUAL INSPECTION TRUKING OF MACHINE							Document No. : SOP/EBI/EN-061-F01A	
								Effective Date :	

Component	Activity							Condition		Remark
	C	I	L	R	S	T	W	B	G	
Pneumatic system										
										
Cleaning condensation water										
Cleaning air filter regulator										
Cleaning solenoid valve										
Air filter in the hose tubing										
Robotic arm										
Lubricant of moving part										
										
Cleaning robotic arm										
Connector cables										
Tighten all bolt and nut										
Motor and gearbox										
Vibration in the robotic arm										
Star wheel										
										
Rubber suction cup										
Fastener starwheel										
L4 Vacuum pump (Becker)										
										
Carbon vanes										
Section gasket										
Section rubber coupling										
Air filter										

 PT ETANA BIOTECHNOLOGIES INDONESIA	FORMULIR / FORM		Document No. : SOP/EBI/EN-061-F01A
	PERAWATAN TERHADAP MESIN AUTOMATIC VISUAL INSPECTION TRUKING		Effective Date :
		MAINTENANCE OF AUTOMATIC VISUAL INSPECTION TRUKING OF MACHINE	

Note :

No	Date	Start	Finish	Duration		Total personnel	Area	Remark
				Hour	Minute			

Remark : C : Cleaning (pembersihan) I : Inspection (pemeriksaan) L : Lubrication (pelumasan) R : Replacement (penggantian) S : Setting (penyetelan) T : Tightened (dikencangkan) W : Welding (pengelasan) G : Good (baik) B : Bad (jelek) * : Critical component	Level : L2 : Every 3 month L4 : Every 1 year	Done by : ()	Checked by : ()
--	--	---	--