 PT ETANA BIOTECHNOLOGIES INDONESIA	STANDARD OPERATING PROCEDURE		Page : 1 of 6
	TITLE : PERAWATAN TERHADAP MESIN BIOSEALER SARTORIUS MAINTENANCE OF BIOSEALER SARTORIUS MACHINE		SOP No. : SOP/EBI/EN-062 Revision : 00 Effective Date : Reviewed Date :
PREPARED BY :		REVIEWED BY :	APPROVED BY :
WENDI RUKMANSYAH (ENGINEERING SPV)		PURNO BUDI KISWANTO (ENGINEERING MANAGER)	VINO SOADUON (PROD MANAGER)
APPROVED BY :		APPROVED BY :	
DEVIKA SRI WULANDARI (QS MANAGER)		SHANDY (QA MANAGER)	

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Provide clear instructions on how to maintenance for Biosealer Sartorius Machine in the production area.

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This SOP is used by the Engineering department as a guideline for carrying out Biosealer Sartorius Machine maintenance in the production area.

3. TANGGUNG JAWAB

RESPONSIBILITY

3.1. Teknisi Engineering *Preventive Maintenance* bertanggung jawab :

Engineering Technician Preventive Maintenance are responsible :

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Performing maintenance duties properly, in accordance with this SOP.

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Perform maintenance according to a predefined schedule.

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Report any damage and repairs that are needed to the Engineering Supervisor.

3.2. Supervisor Engineering bertanggung jawab :

Engineering Supervisor are responsible :

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Make schedule of maintenance periodically.

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Ensure implementation of maintenance is according to the schedule.

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Ensure all measuring instruments involved in the maintenance process are in accordance with the specifications specified (calibrated).

3.2.4. Mengadakan pemeriksaan bulanan terhadap hasil kerja Teknisi Preventive Maintenance.

Perform a monthly review of the workings of the Preventive Maintenance Technicians.

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Provide training to others department which related.

3.3. Manajer Engineering bertanggung jawab :

Engineering Manager is responsible :

3.3.1. Menentukan spare parts pengganti apabila spare parts yang sesuai sudah tidak dijual lagi.

Determine replace parts if suitable spare part are no longer sold.

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Resolve and coordinate necessary actions, if necessary replacement, equipment modification, system changes and prepared change control.

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The Production Machine Operator is responsible for :

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Report if the area is detected by defects / abnormalities that normally impede the operation of the machine.

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Perform daily maintenance such as: tighten ferulle clamps loose bolts, clean the body of the machine / Main Frame at all times so that it does not rust and the surface of the machine remains good, cleaned every machine that is the responsibility, especially on the inside of the crust that is attached.

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Ensure the condition of electrical and electronic equipment that is part of the machine placed in water-prone areas, in order to always avoid splashing water.

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The Calibration Staff of the Quality Assurance Department is responsible for calibrating the components of the measuring instruments used, such as Thermometer, Pressure Gauge, Thermo Controller etc., and conducting regular qualifications.

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Sartorius Biosealer Machine is a tool used in the medical and pharmaceutical industries to protect and maintain the integrity of biomedical products, such as vaccines, medicines, or other biological materials, from contamination or damage. This machine is usually used to seal or close packaging in a very sterile way.

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Autonomous maintenance is an independent engine maintenance carried out by the machine operator. If all this time the operators are only trained to operate the engine, then it's time to be further trained.

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Before implementing preventive maintenance, engineering staff must coordinate with the operator about the area to prepare the area.

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When carrying out maintenance, the related operator must accompany Engineering to avoid operating errors when the machine is running or the operator can undergo autonomous maintenance on the machine.

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Work safety must always be considered during Preventive Maintenance, use PPE (personal protective equipment) according to its designation in conditions during Maintenance. Such as using safety shoes, safety helmets and safety gloves.

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For the replacement of critical components or spare parts that have a direct impact on the product when replaced with different brands and different types but one type, make change control.

6. ALAT DAN BAHAN

TOOLS AND MATERIALS

6.1. ALAT TOOLS

- 6.1.1. 1 set tools box, pada saat melakukan *Preventive Maintenance*.

1 set of toolbox, when performing Preventive Maintenance.

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Level	Component	Activity							Condition		Remark
		C	I	L	R	S	T	W	B	G	
L2											

Socket power											
Button and switch											
Ceramic element heater											
Cleaning cover machine											

Note :

No	Date	Start	Finish	Duration		Total personnel	Area	Remark
				Hour	Minute			

Remark :

C : Cleaning (pembersihan)
 I : Inspection (pemeriksaan)
 L : Lubrication (pelumasan)
 R : Replacement (penggantian)
 S : Setting (penyetelan)
 T : Tightened (dikencangkan)
 W : Welding (pengelasan)
 G : Good (baik)
 B : Bad (jelek)
 * : Critical component

Level :
 L2 : Every 3 month

Done by :

() ()

Checked by :